

# Spectrophotometer Experiment

This document outlines the experiment conducted using a spectrophotometer to measure the absorbance of various materials across different wavelengths. The experiment includes data collection, calibration, and results visualization.

## Experiment Details

1. **Objective**: To measure the absorbance of various materials using a spectrophotometer and to analyze the results.

2. **Materials**:

- Spectrophotometer
- Light Source (LEDs)
- Sample Materials
- Potentiometer
- LDR (Light Dependent Resistor)
- Data Collection Software

3. **Method**:

- Set up the spectrophotometer.
- Calibrate using known reference materials.
- Measure the voltage with and without material at different wavelengths.
- Collect and store data for analysis.

## Results

The results of the experiment will be recorded in a table format, showing the voltage readings with and without the material for various wavelengths.

Additionally, graphs will be generated to visualize the data.

## Conclusion

The experiment successfully demonstrated how absorbance can be measured using a spectrophotometer. The data collected will aid in understanding the absorption characteristics of different materials.